## Alfa Amilase G7

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	CHO
Full Name	Cholesterol
Method	End Point
	Single Reagent
Filter	510
SubFilter	670
Decimal	x.xxx
Unit	mmol/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	1.0
Sample Volume	4
Dilution	1
Normal Low	3.1
Normal High	5.2
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	35 - 36
Linear Range	19.4
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardenha, 110

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## Sinnowa D280

Item	CREA
Full Name	Creatinine
Method	Two Point
	Double Reagent
Filter	510
SubFilter	670
Decimal	xx.xx
Unit	µmol/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.5
Sample Volume	16
Dilution	1
Normal Low	53
Normal High	115
Reagent 1 Volume	320
Reagent 2 Volume	80
Test Point	8 - 15
Linear Range	1330
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	ALP
Full Name	Alk Phos
Method	Kinetic
	Single Reagent
Filter	405
SubFilter	510
Decimal	xxx.x
Unit	U/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.0
Sample Volume	8
Dilution	1
Normal Low	0
Normal High	141
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	5 - 15
Linear Range	700
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	-
Standard Factor	2757**
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

#) Data entry by the user

\*\* ) Factor must be checked by a calibration serum

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	GLU
Full Name	Glucose
Method	End Point Single Reagent
Filter	510
SubFilter	670
Decimal	x.xxx
Unit	mmol/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	1.0
Sample Volume	4
Dilution	1
Normal Low	3.9
Normal High	6.4
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	35 - 36
Linear Range	22.2
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	HDL-C
Full Name	HDL Cholesterol
Method	End Point Double Reagent
Filter	620
SubFilter	670
Decimal	xxx.x
Unit	mg/dL
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	1.0
Sample Volume	5
Dilution	1
Normal Low	35
Normal High	80
Reagent 1 Volume	400
Reagent 2 Volume	130
Test Point	4 - 36
Linear Range	180
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	LDH
Full Name	LDH
Method	Kinetic
	Single Reagent
Filter	340
SubFilter	405
Decimal	xxx.x
Unit	U/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.0
Sample Volume	4
Dilution	1
Normal Low	110
Normal High	248
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	5 - 15
Linear Range	1200
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	-
Standard Factor	16238**
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

#) Data entry by the user

\*\* ) Factor must be checked by a calibration serum

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	LIP
Full Name	Lipase
Method	Kinetic Duoble Reagent
Filter	578
SubFilter	670
Decimal	xxx.x
Unit	U/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.0
Sample Volume	7
Dilution	1
Normal Low	0
Normal High	60
Reagent 1 Volume	320
Reagent 2 Volume	80
Test Point	12 - 19
Linear Range	300
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	MG
Full Name	Magnesium
Method	End Point
	Single Reagent
Filter	546
SubFilter	670
Decimal	x.xxx
Unit	mg/dl
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	1.0
Sample Volume	4
Dilution	1
Normal Low	1.8
Normal High	2.6
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	35 - 36
Linear Range	5.0
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	AST
Full Name	AST
Method	Kinetic
	Single Reagent
Filter	340
SubFilter	405
Decimal	xxx.x
Unit	U/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.0
Sample Volume	40
Dilution	1
Normal Low	0
Normal High	35
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	5 - 15
Linear Range	280
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	-
Standard Factor	1768**
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

#) Data entry by the user

\*\* ) Factor must be checked by a calibration serum

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	ALT
Full Name	ALT
Method	Kinetic
	Single Reagent
Filter	340
SubFilter	405
Decimal	xxx.x
Unit	U/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.0
Sample Volume	40
Dilution	1
Normal Low	0
Normal High	45
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	5 - 15
Linear Range	280
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	-
Standard Factor	1768**
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

#) Data entry by the user

\*\* ) Factor must be checked by a calibration serum

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	TG
Full Name	Triglycerides
Method	End Point Single Reagent
Filter	510
SubFilter	670
Decimal	x.xxx
Unit	mmol/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	1.0
Sample Volume	4
Dilution	1
Normal Low	0.4
Normal High	2.3
Reagent 1 Volume	400
Reagent 2 Volume	-
Test Point	35 - 36
Linear Range	11.3
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

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Rua Cristovão Sardinha, 110

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## Sinnowa D280

Item	UREA
Full Name	UREA
Method	Two Point Double Reagent
Filter	340
SubFilter	405
Decimal	xx.xx
Unit	mmol/L
Blank Medium	Reagent
Blank Value	-
Blank Low	0
Blank High	2.0
Sample Volume	3
Dilution	1
Normal Low	2.8
Normal High	7.2
Reagent 1 Volume	320
Reagent 2 Volume	80
Test Point	8 - 12
Linear Range	50.0
Assistant Start	NONE
Assistant End	NONE
Range Low	-
Range High	-
Number of Standards	1
Standard Factor	-
Correct Factor	-
Standard Concentration	#
Sample Storage Time (2 – 8 C)	7 days

#) Data entry by the user