**VITROS 5,1 FS**

TEST	CRP
TESTTYP	end-point
SAMPLE TYPE	serum
CALIBRATION TYPE	Logit / Log4
TEMPLATE	EPT1 R1 – S – R2
CAL. No.	6
DIL. FACTOR STANDARDI	1.0
DILUTION	saline
UNIT	mg/dl
SLOPE	1.00
Y_AXIS INTERCEPT	0.00
REFERENCE RANGE	#
BLANK ABS	-0.200 – 2.700
MONOTONIC	increase
REACTION RANGE	-3.000 – 3.000
CAL.CURVE ADJUSTM:	0.99
MIN / MAX VALUE	0.01 / 30.0
CALIBRATION CONC:	*.0 - # - # - # - # - #
<b>PROTOCOL STEPS</b>	
1. REAGENT VOL.(µL)	200
2. INCUBATION (sec)	0.00
3. SAMPLE (µL)	13
4. INCUBATION (sec)	60
5. MEASUREMENT wavelength (nm)	340
6. INCUBATION (sec)	4.75
7. REAGENT VOL. (µL)	40
8. INCUBATION (sec)	313.50
9 MEASUREMENT wavelength (nm)	340

#) Data entry by the user

\*) Use 0.9% NaCl as zero calibrator

**Order information**

**Cat. No. 1 7002 . . . . .**

**Notes**

1. Please refer to the package insert for CRP FS for detailed information about the test on the following:

Clinical Relevance  
Method and Principle  
Composition and Stability of the Reagents  
Specimens  
Calibrators and Controls  
Performance Characteristics regarding  
- Measuring Range  
- Specificity/Interferences  
- Sensitivity/Limit of Detection  
- Precision (Reproducibility, Repeatability)  
- Method Comparison  
Reference Ranges  
Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
3. Manufactured by  
DiaSys Diagnostic Systems GmbH  
Alte Strasse 9, 65558 Holzheim, Germany

**\*\*This application proposal is for guidelines only. To avoid misinterpretation measured results have to be validated and assessed with caution.**

# APPLICATION NOTE

## VITROS 5.1 FS

### C3C KIT

#### 1. Reagent Preparation

Sample : Ready for use.  
 Reagent 1 : Buffer, ready for use.  
 Reagent 2 : Latex, ready for use.  
 Calibration : Dilute Protein Standard High successively 1:2 in saline 9 g/L in order to set up a calibration curve. Alternatively use the ready for use Protein Standard Set. Use 9 g/L NaCl as zero point.

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\*\* see insert

### 3. Ordering information

107C004 1 x 10 mL Antiserum  
2 x 25 mL Buffer  
MPS/STH Protein Standard High, 5 mL  
MPS/STS-5X1 Protein Standard Set, 5x1 mL  
139F003 Immunology Control Low, 1mL  
139F002 Immunology Control High, 1mL

# APPLICATION NOTE

## VITROS 5.1 FS

### C4C KIT

#### 1. Reagent Preparation

Sample : Ready for use.  
 Reagent 1 : Buffer, ready for use.  
 Reagent 2 : Antiserum, ready for use.  
 Calibration : Dilute Protein Standard High successively 1:2 in saline 9 g/L in order to set up a calibration curve. Alternatively use the ready for use Protein Standard Set. Use 9 g/L NaCl as zero point.

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\*\* see insert

### 3. Ordering information

108C004      1 x 10 mL Antiserum  
                   2 x 25 mL Buffer  
 MPS/STH     Protein Standard High, 5 mL  
 MPS/STS-5X1 Protein Standard Set, 5x1 mL  
 139F003     Immunology Control Low, 1mL  
 139F002     Immunology Control High, 1mL

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