## Ácido Úrico

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Kovalent do Brasil Ltda.

Rua Cristovão Sardenha, 110

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Test	Uric Acid
Units	mg/dL
Analysis Mode	Final Point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Bichromatic
Primary Wavelength	505 nm
Reference Wavelength	670 nm
Volume of Sample	8 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	600 s
Incubation Time 2	0
Stabilisation Time	15 s
OPTIONS	
Linearity Limit	20 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	8.20
Lower Reference Range	2.30
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Albumin
Units	g/dL
Analysis Mode	Final point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	670 nm
Reference Wavelength	-
Volume of Sample	5 µL
Volume of Reagent 1	500 µL
Volume of Reagent 2	0 µL
Incubation time 1	300 s
Incubation Time 2	0
Stabilisation Time	5 s
OPTIONS	
Linearity Limit	6 g/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	5.2
Lower Reference Range	3.5
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Amylase
Units	U/L
Analysis Mode	Kinetic
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Value of the factor	4554
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	405 nm
Reference Wavelength	-
Volume of Sample	8 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	50 s
Incubation Time 2	60 s
Stabilisation Time	30 s
OPTIONS	
Linearity Limit	2000 U/L
Blank Absorbance Limit	1 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	100 U/L
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Test	TDil
Units	mg/dL
Analysis Mode	Differential
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	546 nm
Reference Wavelength	-
Volume of Sample	40 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	100 µL
Incubation time 1	180 s
Incubation Time 2	300 s
Stabilisation Time	10 s
OPTIONS	
Linearity Limit	10 mg/dL
Blank Absorbance Limit	1 A
Kinetic Blank Limit	.2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	0.2 mg/dL
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	TBil
Units	mg/dL
Analysis Mode	Differential
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	546 nm
Reference Wavelength	-
Volume of Sample	10 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	100 µL
Incubation time 1	300 s
Incubation Time 2	300 s
Stabilisation Time	10 s
OPTIONS	
Linearity Limit	30 mg/dL
Blank Absorbance Limit	1 A
Kinetic Blank Limit	.2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	1.2 mg/dL
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Calcium
Units	mg/dL
Analysis Mode	Final Point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	670 nm
Reference Wavelength	nm
Volume of Sample	4 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	300 s
Incubation Time 2	0
Stabilisation Time	5 s
OPTIONS	
Linearity Limit	25 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	10.3 mg/dL
Lower Reference Range	8,6 mg/dL
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	CK MB
Units	U/L
Analysis Mode	Kinetic
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Value of the factor	8254
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	0
Reading	Monochromatic
Primary Wavelength	340 nm
Reference Wavelength	-
Volume of Sample	16 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	300 s
Incubation Time 2	120 s
Stabilisation Time	30 s
OPTIONS	
Linearity Limit	2000 U/L
Blank Absorbance Limit	0.5 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	25 U/L
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	CK NAC
Units	U/L
Analysis Mode	Kinetic
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Value of the factor	4127
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	0
Reading	Monochromatic
Primary Wavelength	340 nm
Reference Wavelength	-
Volume of Sample	16 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	120 s
Incubation Time 2	60 s
Stabilisation Time	10 s
OPTIONS	
Linearity Limit	2000 U/L
Blank Absorbance Limit	0.5 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	171 U/L
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Cholesterol
Units	mg/dL
Analysis Mode	Final Point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Bichromatic
Primary Wavelength	505 nm
Reference Wavelength	700 nm
Volume of Sample	4 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	600 s
Incubation Time 2	0
Stabilisation Time	10 s
OPTIONS	
Linearity Limit	750 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	200 mg/dL
Lower Reference Range	# mg/dL
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Test	Iron
Units	µg/dL
Analysis Mode	Dif
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	578 nm
Reference Wavelength	-
Volume of Sample	50 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	100 µL
Incubation time 1	300 s
Incubation Time 2	300s
Stabilisation Time	5 s
OPTIONS	
Linearity Limit	1000 µg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	5.2
Lower Reference Range	3.5
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Glucose
Units	mg/dL
Analysis Mode	Final Point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Bichromatic
Primary Wavelength	505 nm
Reference Wavelength	670 nm
Volume of Sample	4 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	600 s
Incubation Time 2	0
Stabilisation Time	5 s
OPTIONS	
Linearity Limit	400 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	115 mg/dL
Lower Reference Range	75 mg/dL
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	500 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	HDL-Col
Units	mg/dL
Analysis Mode	Final Point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	0
Reading	Bichromatic
Primary Wavelength	505 nm
Reference Wavelength	700 nm
Volume of Sample	4 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	100 µL
Incubation time 1	300 s
Incubation Time 2	300 s
Stabilisation Time	10 s
OPTIONS	
Linearity Limit	180 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	35 mg/dL
Lower Reference Range	# mg/dL
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	AST/GOT
Units	U/L
Analysis Mode	Kinetic
Associated constituents	#
Reaction Type	Decrease
Sample replicates	1
Value of the factor	1745
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	340 nm
Reference Wavelength	-
Volume of Sample	40 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	50 s
Incubation Time 2	60 s
Stabilisation Time	30 s
OPTIONS	
Linearity Limit	600 U/L
Blank Absorbance Limit	1 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	35 U/L
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Test	ALT/GPT
Units	U/L
Analysis Mode	Kinetic
Associated constituents	#
Reaction Type	Decrease
Sample replicates	1
Value of the factor	1745
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Monochromatic
Primary Wavelength	340 nm
Reference Wavelength	-
Volume of Sample	40 µL
Volume of Reagent 1	400 µL
Volume of Reagent 2	0 µL
Incubation time 1	50 s
Incubation Time 2	60 s
Stabilisation Time	30 s
OPTIONS	
Linearity Limit	600 U/L
Blank Absorbance Limit	1 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	41 U/L
Lower Reference Range	#
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	400 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Triglycerides
Units	mg/dL
Analysis Mode	Final Point
Associated constituents	#
Reaction Type	Increase
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	2
Reading	Bichromatic
Primary Wavelength	546 nm
Reference Wavelength	670 nm
Volume of Sample	5 µL
Volume of Reagent 1	500 µL
Volume of Reagent 2	0 µL
Incubation time 1	600 s
Incubation Time 2	0
Stabilisation Time	10 s
OPTIONS	
Linearity Limit	1000 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	200 mg/dL
Lower Reference Range	# mg/dL
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	500 µL

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Kovalent do Brasil Ltda.

Rua Cristovão Sardinha, 110

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Test	Urea
Units	mg/dL
Analysis Mode	Fixed Time
Associated constituents	#
Reaction Type	Decrease
Sample replicates	1
Calibrator	#
Rep Calibrator and Blank	2
Number of Calibrators	#
Calibrator 1	*
Decimals	1
Reading	Monochromatic
Primary Wavelength	340 nm
Reference Wavelength	-
Volume of Sample	5 µL
Volume of Reagent 1	500 µL
Volume of Reagent 2	0 µL
Incubation time 1	50 s
Incubation Time 2	60
Stabilisation Time	30 s
OPTIONS	
Linearity Limit	300 mg/dL
Blank Absorbance Limit	0 A
Kinetic Blank Limit	2 A/min
High Factor Limit	#
Low Factor Limit	#
Upper Reference Range	43
Lower Reference Range	17
Number of Controls	#
Control Types	#
Control Replicates	#
Washing Volume	500 µL